







www.CPAPDOCTOR.ORG
FOR PATIENTS AND CLINICIANS

ResMed AirFit F30 Full Face Mask

Cleaning/Replacement Guidelines

	ResMed Product#	HCPSC Code	Replacement Schedule
 AirFit F30 mask (with headgear)	AirFit F30 mask with headgear, with cushion size: <ul style="list-style-type: none"> • S - 64100 • M - 64101 	A7034 + A7035 (Mask + Headgear)	
 AirFit F30 Frame System	AirFit F30 frame with cushion size: <ul style="list-style-type: none"> • S - 64155 • M - 64156 	A7034 (Mask)	Replace every 3 months.
 AirFit F30 Headgear	AirFit F30 Headgear: <ul style="list-style-type: none"> • 64161 	A7035 (Headgear)	Replace every 6 months.
 AirFit F30 Cushion	AirFit F30 cushion size: <ul style="list-style-type: none"> • S - 64150 • M - 64151 	A7033 (Cushion)	Replace at least every 3 months. May replace every month if needed.
Cleaning Schedule	Wash weekly with warm water and mild soap and air dry. May rinse frame and cushions daily if needed.		

DME Provider:
Phone Number: