







www.CPAPDOCTOR.ORG
FOR PATIENTS AND CLINICIANS

ResMed AirCurve 10

Cleaning & Replacement Guidelines

	ResMed Product #	HCPCS Code	Cleaning Schedule	Replacement Schedule
 AirCurve 10		E04070 + E0562 (BiLevel without backup rate) or E04071 + E0562 (BiLevel with backup rate)	Wipe weekly the exterior surface with a dry cloth.	Every <u>5 years</u> .
 Standard Air Filter or Hypoallergenic Air Filter*	36851 (2-pack) 36857 (2-pack)	A7038 (Disposable filter)	Check weekly for dirt or blockage. Do NOT wash.	Replace at least every <u>6 months</u> . May be replaced every 2 weeks if needed.
 SlimLine Tubing or ClimateLineAir Heated Tube	36810 37296	A7037 (Standard tubing) A4604 (Heated tubing)	Clean weekly with warm water and mild soap and air dry.	Replace every <u>3 months</u> .
 HumidAir Water Tub	37299	A7046	Clean weekly with warm water and mild soap and air dry.	Replace every <u>6 months</u> .

Standard air filter efficiency > 75% for ~7 micron dust.

*Hypoallergenic air filter > 98% for ~7-8 micron, > 80% for 0.5 micron dust.

DME provider:
Phone number: