



www.CPAPDOCTOR.ORG
FOR PATIENTS AND CLINICIANS

Breas Z2 Auto

Cleaning & Replacement Guidelines

	Breas Product #	HCPCS Code	Cleaning Schedule	Replacement Schedule
 Z2 Auto		E0601 (CPAP)	Wipe weekly the exterior surface with a damp cloth.	Every 5 years .
 Z1/Z2 Polyester Filter	Z1/Z2 Filter (2-pack): • 005788	A7038 (Disposable filter)	Check weekly for dirt or blockage. Do NOT wash.	Replace at least every month . May be replaced every 2 weeks if needed.
 Slim Tube, 4 ft.	Slim Tube (4 ft.): • 005797	A7037 (Standard tubing)	Clean weekly with warm water and mild soap and air dry.	Replace every 3 months .
 Heat Moisture Exchange (HME)	Heat Moisture Exchange (HME): • HD60-2010			Replace every week .

DME Provider:
Phone Number: