







www.CPAPDOCTOR.ORG
FOR PATIENTS AND CLINICIANS

Philips Amara View Full Face Mask

Cleaning/Replacement Guidelines

	Philips Product#	HCPCS Code	Replacement Schedule
 <p>Amara View Full Face Mask (with headgear)</p>	<p>Amara View mask with headgear, with cushion size:</p> <ul style="list-style-type: none"> • Fitpack (S, M & L) - 1090670 • S - 1090622 • M - 1090623 • L - 1090624 	A7034 + A7035 (Mask + Headgear)	
 <p>Amara View Full Face Mask</p>	<p>Amara View mask with cushion size:</p> <ul style="list-style-type: none"> • S - 1090632 • M - 1090633 • L - 1090634 	A7034 (Mask)	Replace every 3 months.
 <p>Amara View Headgear</p>	<p>Amara View Headgear size:</p> <ul style="list-style-type: none"> • Standard - 1090697 • Large - 1090696 	A7035 (Headgear)	Replace every 6 months.
 <p>Amara View Cushion</p>	<p>Amara View Cushion size:</p> <ul style="list-style-type: none"> • S - 1090692 • M - 1090693 • L - 1090694 	A7033 (Cushion)	Replace at least every 3 months. May replace every month if needed.
Cleaning Schedule	<p>Wash weekly with warm water and mild soap and air dry. May wash non-fabric parts daily if needed.</p>		

DME Provider:
Phone Number: