



www.CPAPDOCTOR.ORG
FOR PATIENTS AND CLINICIANS

Fisher & Paykel Vitera Full Face Mask

Cleaning/Replacement Guidelines

	Fisher & Paykel Product#	HCPCS Code	Replacement Schedule
 Vitera Full Face Mask (with headgear)	Vitera full face mask with headgear, with cushion size: • Fit Pack (S, M & L) - VIT1SMLA • S - VIT1SA • M - VIT1MA • L - VIT1LA	A7034 + A7035 (Mask + Headgear)	
 Vitera Full Face Mask	Vitera full face mask with cushion size: • S - 400VIT131 • M - 400VIT132 • L - 400VIT133	A7034 (Mask)	Replace every 3 months .
 Vitera Headgear	Vitera Headgear size: • S - 400VIT121 • M/L - 400VIT122	A7035 (Headgear)	Replace every 6 months .
 Vitera Seal Cushion	Vitera Seal Cushion size: • S - 400VIT114 • M - 400VIT115 • L - 400VIT116	A7033 (Cushion)	Replace at least every 3 months. May replace every month needed.
Cleaning Schedule	Wash weekly with warm water and mild soap and air dry. May wash mask excluding the headgear daily if needed.		

DME Provider:
Phone Number: