







www.CPAPDOCTOR.ORG  
FOR PATIENTS AND CLINICIANS

## ResMed AirFit F30i Full Face Mask

### Cleaning/Replacement Guidelines

	ResMed Product#	HCPCS Code	Replacement Schedule
 <b>AirFit F30i mask</b> (with headgear)	AirFit F30i mask, <b>standard</b> frame with cushion size: <ul style="list-style-type: none"> <li>• S - <b>63331</b></li> <li>• M - <b>63332</b></li> <li>• W - <b>63333</b></li> </ul> AirFit F30i mask, <b>small</b> frame with cushion size: <ul style="list-style-type: none"> <li>• S - <b>63330</b></li> </ul> AirFit F30i mask, <b>large</b> frame with cushion size: <ul style="list-style-type: none"> <li>• M - <b>63334</b></li> </ul>	<b>A7034 + A7035</b> (Mask + Headgear)	
 <b>AirFit F30i</b> <b>Frame System</b>	AirFit F30i mask, <b>standard</b> frame with cushion size: S - <b>63359</b> , M - <b>63360</b> SW - <b>63361</b> , W - <b>63362</b> AirFit F30i mask, <b>small</b> frame with cushion size: S - <b>63355</b> , M - <b>63356</b> SW - <b>63357</b> , W - <b>63358</b> AirFit F30i large frame with cushion size: S - <b>63363</b> , M - <b>63364</b> SW - <b>63365</b> , W - <b>63366</b>	<b>A7034</b> (Mask)	Replace every <b>3 months.</b>
 <b>AirFit F30i</b> <b>Headgear</b>	AirFit F30i Headgear: <ul style="list-style-type: none"> <li>• <b>63372</b></li> </ul>	<b>A7035</b> (Headgear)	Replace every <b>6 months.</b>
 <b>AirFit F30i</b> <b>Cushion</b>	AirFit F30i cushion size: <ul style="list-style-type: none"> <li>• S - <b>63350</b></li> <li>• M - <b>63351</b></li> <li>• SW - <b>63352</b></li> <li>• W - <b>63353</b></li> </ul>	<b>A7033</b> (Cushion)	Replace at least every <b>3 months.</b>  May replace every month if needed.
<b>Cleaning Schedule</b>	Wash <b>weekly</b> with warm water and mild soap and air dry. May rinse cushions daily if needed.		

DME Provider:  
Phone Number: